

Options for NLP 2000 System

Localized Environmental Chamber (EC-0110-00)

User controlled temperature and humidity within the printing environment.

- Software controlled temperature & humidity
- Heating temperature range: (ambient + 20°C)
- Cooling temperature range: (ambient - 2°C)
- Temperature stability: +/- 0.5 °C
- Humidity range: (10-90% RH)
- Humidity stability: +/- 5% RH



2D nano PrintArray Option (NLP-0100-01)

The 2D nano PrintArray™ Chip provides a high-throughput solution to flexibly pattern nanoscale features. The 2D nano PrintArray retains the direct write, high resolution, ambient deposition, and chemical and material flexible attributes of DPN, while multiplying the desired pattern 55,000 times across a 1 cm² area.


Extended Service Package (NLP-0310-01)

Receive an additional 1-year of service and warranty beyond the standard 1-year.

Integrated Vibration Isolation Option (ATB-0100-00)

The NLP 2000 System has an integrated passive vibration isolation system.

- Vistek bearings isolate horizontal and vertical vibration components
- No external air source required
- System footprint: 30" (76 cm) L x 22" (56 cm) W



Compatible with the NLP 2000 System

Learn more about NanoInk products and services at www.nanoink.net. Or call us at 847-679-NANO (6266).